

Mechanical profilometer (stylus based)

*Refer to the images attached for actual position, placement and settings.

1. Power on the computer (under the table)
2. Logon info is on computer monitor frame
3. Turn on machine (behind computer monitor)
4. Open 'zworld' application. As long as "C" is being transmitted, the connection is good and the data collection software can be opened.
5. Open 'Dektak' software. (not Dektak offline)
6. Click on tower down stylus up (found on the top left bar with little images, use cursor to move over them and find tower down -stylus up). Make absolutely sure that the stylus is above the sample and not the stage.
7. Change illumination levels to see the sample.
8. Calibrate zoom and focus of the microscope attached above the stylus.
9. The stage can be rotated.
10. The scan axis is in and out or vertically displayed on the screen.
11. Go to window-scan routine.
12. Modify scan time, scan length, tip pressure
13. Tip pressure you can start at 6mg
14. Click run currently active scan routine.
15. For a map scan, go to automation program
16. Select a file to save the data to
17. Click run automation program
18. Stage will retract automatically when closed.
19. Use 'vision' software to export dektak data to a useable format.
20. Turn off by closing 'Dektak' application, turn off tool (switch off), then shutdown computer.

