

E-Beam Lithography

Friday, April 4, 2025 1:39 PM

<https://cns1.rc.fas.harvard.edu/facilities/docs/2010-07-09%20EBL%20talk.pdf>

- Steps:
- GDS Marks
 - Beamer
 - WeCas
 - SEM
 - WeCas

GDS Marks → Beamer

- 1) Get your GDS file with marks (layers)
- 2) Edit your GDS in Layout editor mode to get the desired layer for writing.
- 3) Save it

4) open Beamer App

- 5) Create Process flow / Load existing flow → For 10T's
 • Load GDS
 • Heating
 • Transform
 • PEC ⇒ use right stimulation
 • Select dose rate in the last process
 • Save .con file

6) Go to WEcas in laptop


- Edit Schedule
- Load .con file
- Check estimated time
- Dose rate calculation
- Get estimate and Book tool

7) Get the Spin coated Sample

8) Load the Sample

- Check exchange position in SEM app
- On the machine check if Sample Stage in "ic"
- Open and Load the Sample
- Close and Press "Evac". and wait till pressure Vacuum
- open gate
- Push rot in and unlock the stage and pull back
- Close gate

9) SEM app

-  click to Recall Beam
- Richfield 3nA 120 c/s
- Isolation Valve Button click
- Faraday Cup
- Zoom out (Max) and center the green dots to center the circle
- Press Measure Current
- For 3nA if it not around 3, Click adjust to increase or decrease the current.
- Then goto Reference position
- Find a dust and focus
- Move stage often to avoid charging
- Find gold particles
- Focus and sigma  alternatively
- After seeing the particles clear
- Move Faraday Cup and Measure current / adjust
- Then check Reference again for focus and clarity

10) WEcas app

- Edit Schedule
- Load .con file and next
-  Move pattern
- Measure Sample inclination and next (the laser displacement meter should be 0.000)
- Field Correction (OK)
- Calc.
- Exposure
- Wait till a few field is written to check for inclination error.

Move stage to left bottom corner before measuring inclination

Look at the z preset value (maybe change it if you get the laser displacement meter not being in 0.0000)

*** DOSE CALCULATION ***

DOT MAP	60000	dots
AREA	150	micron ²
REG DOSE	650.000	micro C/cm ²
BEAM CURRENT	2.00e-010	A
DOSE TIME	0.2031	micro sec.
CALC		EXIT





Flow chart of Elionix operating procedure

BEAM CURRENT	12.00e-010	A
DOSE TIME	0.2031	micro sec.
CALC		EXIT


https://www.seas.upenn.edu/~nanosop/Elionix_Alignment_SOP.htm

6. Confirm the coordinates of Reg-A and Reg-B marks using SEM

1. Open the "Stage Memory" window by clicking the "Stage Memory" icon on the Tool bar.
2. Type the estimated coordinate of Reg-A mark in the "Stage Memory" window.
3. Click the "Move" button.



4. After the stage moves, open the beam shutter.
5. Find the Reg-A mark on the SEM image.
6. Check the "Cross Mark" in the "Scan" window, if not.
7. Move the center of the Reg-A mark to the center of the cross mark, using the mouse (select "Move Stage (Center)" from the drop down menu).
8. Increase the magnification (x2000 - x100000), and then move the center of the Reg-A mark to the center of the cross mark.
9. Record the coordinate of the Reg-A mark shown on the "Position" window.
 - This is the accurate coordinate of the Reg-A mark.
 - Ensure that the slow mode in the "Scan" window is selected.



10. Close the beam shutter.
11. Move the stage to the estimated coordinate of Reg-B mark in the same manner as above.
12. After the stage moves, open the beam shutter.
13. Move the center of the Reg-B mark to the center of the cross mark, and record the coordinate of the Reg-B mark in the same manner as above.
 - This is the accurate coordinate of the Reg-B mark.
14. Close the beam shutter.

IMPORTANT:

- For example, it is assumed that the designed coordinates of the center, Reg-A, and Reg-B are (0,0), (0,-2), and (0,2) in mm, respectively. The sample size is 15 mm x 15 mm. Then, the absolute coordinates of the center, Reg-A, and Reg-B on the tool are (84.3, 85.5), (84.3, 83.5), and (84.3, 87.5) in mm, respectively. The relative relationship must be precise. The Elionix will adjust the real coordinates during exposure procedure, as shown later.




<https://www.yumpu.com/en/document/read/19278682/exposure-procedure/18>

Place Fields(cc7/cb7 file name)

ELIONIX
Confidential

◆ Manual Placement, OR Auto Placement

- Auto Placement Need...
- Field Name (ccc/cbc file name)
- Manual Placement Need...
- Field Position (X,Y)
- Field Name (ccc/cbc file name)



2013/2/28

Page 12

Save CON file (.con) and next step

ELIONIX
Confidential

- ◆ Save Condition file (ccc/cbc also saved)
- ◆ Move onto Schedule file preparation



2013/2/28

Page 13